




FILTER TECHNOLOGY

FLAMETIP

Pipette Micro Filter Tips



Description

Natural, with filter, sterilized in rack

Materials

Tips: PP (Polypropylene)

Color: Natural

Rack box: PP (Polypropylene)

Purpose

Use for liquid transfer, zero retention filter.

Features

- ◆ FLAMETIP are manufactured from super clear high quality Polypropylene
- ◆ The surfaces of the tips are produced through a special process, this process makes the tip inner surface become super hydrophobic
- ◆ Preferred accessories for most brand micropipettor
- ◆ Tips with PP filter are also available
- ◆ Packaged in re-sealable plastic bags or extra-rigid autoclavable racks
- ◆ Every rack or case is printed with lot No. for quality traceability
- ◆ Non-Pyrogenic and validated per FDA guidelines on LAL testing for medical devices and company guidelines. The acceptance level for product is less than 0.5 EU/ml
- ◆ DNase/RNase-free
- ◆ Sterilized by gamma irradiation SAL 10-6(1S011137)
- ◆ Shelf life: 3 years after month of production
- ◆ Manufactured in a Class 100,000 cleanroom environment
- ◆ Manufactured under ISO13485:2016 and ISO9001:2015 quality management system
- ◆ This product has been tested and is free of any detectable Nucleic Acid.
- ◆ All material are certified metal-free which compliance with US Pharmacopoeia guidelines.
- ◆ This product has been tested and is free of ATP

FLAMETIP Pipette Micro Filter Tips Ordering informations

CODE	DESCRIPTION	STERILE	QTY/PACK	QTY/BOX	QTY/CASE
FTS252010R1	10 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS233010R1	10 µl , long, natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252020R1	20 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252100R1	100 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS231200R1	200 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS233300R1	300 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252000R1	1000 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS371000R1	1000 µl , long, natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920