

Bacteriological Petri Dishes



FILTER TECHNOLOGY

Petri Dishes, Plastic Material

Single-use Petri dishes for microbiological procedures, as well suspension culture or plant cell culture

- Non-cytotoxic virgin polystyrene, excellent optical clarity
- Engineered for optimum flatness to provide uniform agar thickness
- Single or more vents are available: no-vent for anaerobic and long-term work, single vent for limiting gas exchange, multi vents for improving gas exchange



Compartmented Dishes



Square Dishes



Contact Plates

General petri Dishes

Φ90 and 150mm Petri dishes are manufactured in accordance with standard ISO 24998:2008, ideal for use with most automatic filling machine

P/N	DESCRIPTION	PACKING
PETRIDPC035151NVS	Φ35×15mm, 1Room	10pcs/polybag, 300pcs/carton
PETRIDPC049131NVS	Φ49×13mm, 1Room	10pcs/polybag, 1,000pcs/carton
PETRIDPC055151NVS	Φ55×15mm, 1Room	10pcs/polybag, 1,000pcs/carton
PETRIDPC0601513VS	Φ60×15mm, 1Room/3 vents	10pcs/polybag, 1,000pcs/carton
PETRIDPC070151NVS	Φ70×15mm, 1Room	20pcs/polybag, 800pcs/carton
PETRIDPC0901513VS	Φ90×15mm, 1Room/3 vents	10pcs/polybag, 500pcs/carton
PETRIDPC1501513VS	Φ150×15mm, 1Room/3 vents	10pcs/polybag, 200pcs/carton

Compartmented Dishes

Ideal for different media on same dish, or meet the requirement of saving incubator space

P/N	DESCRIPTION	PACKING
PETRIDPC0901523VS	Φ90×15mm, Φ2Rooms/3 vents	10pcs/polybag, 500pcs/carton
PETRIDPC0901533VS	Φ90×15mm, Φ3Rooms/3 vents	10pcs/polybag, 500pcs/carton

Square Dishes

Noncompartmented, ideal for antibiotic sensitivity testing

P/N	DESCRIPTION	PACKING
PETRIDPS130151NVS	130×130×15mm, 1Room	1pcs/polybag, 150pcs/carton

Contact Plates

10×10mm counting net, alphanumeric, facilitates colony counting and locating

P/N	DESCRIPTION	PACKING
PETRIDPC0651513VS	Φ65×15mm, 1Room/3 vents	10pcs/polybag, 1,000pcs/carton

Petri Dishes, Glass Material

3.3 Borosilicate glass, withstand repeated sterilization (wet or dry)



P/N	DESCRIPTION	PACKING
PETRIDGC060151NVS	Ø60mm x15mm	288pcs/carton
PETRIDGC075151NVS	Ø75mm x15mm	288pcs/carton
PETRIDGC090181NVS	Ø90mm x18mm	192pcs/carton
PETRIDGC100201NVS	Ø100mm x20mm	192pcs/carton

Petri Dish Carrier Rack

Designed for holding and carrying plastic or glass petri dishes

- Made of STAINLESS STEEL Material
- Autoclavable and resistant to universal chemicals
- Highly polished surface, no edges and sharps
- Lightweight, easy to load, safe for transporting

Packing: 1pc/inner box, 5pcs/carton

P/N	DIM.(mm)	DESCRIPTION
PETRIDC340110380330	340×110×380	3 sections, holds 30 petri dishes

