



LIFE SCIENCE PRODUCT

PLASTIC CONSUMABLES
RNase-free/DNase-free/Non-pyrogenic





INDEX

| | Pag. |
|--------------------------------------|-------|
| Cryogenic Tube | 1-3 |
| Screw Cap Micro Tubes | 4-5 |
| Transport Tube with Swab Capture Cap | 6-7 |
| Pre-Filling Transport Tube | 8-9 |
| Pipette Micro Filter Tips | 10-21 |






Cryogenic Tube

Cryogenic Tube

Do not use Cryogenic tubes in the liquid phase of liquid nitrogen. These cryogenic tubes are suited for research and general laboratory storage of cell cultures, bacteria and other samples in the vapor phase of liquid nitrogen.


- Produced from virgin PP material, the cap with a silicone gasket mounted in cap for additional leak-proof assurance
- Large writing areas for more sample information printed graduation for accurate measurement
- Sterilized and certified non-pyrogenic and RNase-free and DNase-free
- Self standing, star shaped and round bottom are available



External Thread (1ml/2ml/4ml/5ml)

External Thread (2ml/5ml)

External Thread (2ml/5ml)



Cryogenic Tube, External Thread, PP Material

External thread type can minimize liquid retention.

| P/N | TYPE | VOL.(mL) | DIM.(H X OD)(mm) | GRADUATION(mL) | PACKING |
|----------------|------|----------|------------------|---------------------|--|
| PLAC4004X1749S | A | 1.0 | 41 x12.4 | 0.5/1.0 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1725S | A | 2.0 | 48.5 x12.4 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1744S | B | 2.0 | 48.5 x12.4 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1745S | C | 2.0 | 48.5 x12.4 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1754S | A | 4.0 | 70.3 x12.4 | 1.0/2.0/3.0/3.6 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1726S | A | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1746S | B | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1737S | C | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |

Cryogenic Tube, Internal Thread, PP Material

Internal thread type can minimize liquid retention and suit for storage in the vapor phase of liquid nitrogen.

| P/N | TYPE | VOL.(mL) | DIM.(H X OD)(mm) | GRADUATION(mL) | PACKING |
|----------------|------|----------|------------------|---------------------|--|
| PLAC4004X1721S | A | 1.0 | 41 x12.3 | 0.5/1.0 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1722S | A | 2.0 | 49 x12.3 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1741S | B | 2.0 | 49 x12.3 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1742S | C | 2.0 | 49 x12.3 | 0.5/1.0/1.8 | 50pcs/zip-lock bag, 500pcs/inner box, 2000pcs/carton |
| PLAC4004X1723S | A | 4.0 | 70.3 x12.3 | 1.0/2.0/3.0/3.6 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1738S | A | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1739S | B | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |
| PLAC4004X1740S | C | 5.0 | 91.3 x12.4 | 1.0/2.0/3.0/4.0/4.5 | 25pcs/zip-lock bag, 250pcs/inner box, 1000pcs/carton |

Screw Cap Micro Tubes

Screw Cap Micro Tubes

The secure screw cap with O-ring seal in conjunction with the robust design makes the screw cap micro tube indispensable in your lab. They are ideal for centrifugation, transport, storage, and short-term cold storage.

- Made from 100% virgin polypropylene to help ensure maximum clarity for biological-based work
- Uniform wall thickness to withstand RCF 12000 xg
- Knurled Cap helps you achieve a firm grip when opening and closing
- Knurls lock tube into specialty racks for single-handed operation
- Screw cap micro tubes are available with the option of multicolored caps, attached or unattached cap

Packing: 100pcs/bag, 1000pcs/carton

| P/N | VOL.(ML) | STERILE |
|----------------|----------|---------|
| PLAC4004X1835S | 1.5 | Sterile |
| PLAC4004X1836S | 2 | Sterile |



Transport Tube with Swab Capture Cap


Transport Tube with Swab Capture Cap



- Tube in polypropylene /Cap in polyethylene
- Transport tubes are designed for the safe storage and transport of biological materials.
- The polyethylene cap feature a swab capture design to ensure the swab can be hold securely inside the cap, making it ideal for specimen collection and protection from contamination.
- The transparent polypropylene tubes allow unobstructed viewing of the sample.
- The tubes are shatter-proof and may be exposed to freezing temperatures.
- The tubes feature a frosted writing surface to make sample identification more convenient.
- Offered in five popular sizes




Swab Capture Design



| | P/N | VOLUME | DIM(MM) | GRADUATIONS | TUBE/CAP | PACKING |
|---|-----------------|--------|-----------|-------------------|-----------|---|
| 1 | PLAC4004X1514NS | 5ml | 16 x 56.5 | Molded graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1515NS | 5ml | 16 x 56.5 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 2 | PLAC4004X1562NS | 5ml | 16 x 60 | Molded graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1563NS | 5ml | 16 x 60 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 3 | PLAC4004X1512NS | 8ml | 16 x 85 | Molded graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1513NS | 8ml | 16 x 85 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 4 | PLAC4004X1579NS | 10ml | 16 x 94 | Molded graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1580NS | 10ml | 16 x 94 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 5 | PLAC4004X1590NS | 10ml | 16 x 102 | Non graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1592NS | 10ml | 16 x 102 | Non graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 6 | PLAC4004X1597NS | 7ml | 13 x 80 | Non graduation | Assembled | 100pcs/bag, 1000pcs/carton |
| | PLAC4004X1596NS | 7ml | 13 x 80 | Non graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |

*Dimensions for tube only



| | A(mm) | B(mm) | C(mm) | D(mm) |
|---|-------|-------|-------|-------|
| 1 | 21 | 16 | 58 | 56.5 |
| 2 | 21 | 16 | 62 | 60 |
| 3 | 21 | 16 | 89.5 | 85 |
| 4 | 21 | 16 | 96 | 94 |
| 5 | 21 | 16 | 102 | 100 |
| 6 | 17 | 13 | 80 | 82 |

Pre-Filling Transport Tube

Pre-Filling Transport Tube



- Tube in polypropylene /Cap in polyethylene
- The standard pre-filling transport tube with the perfect seal is designed for pre-filling of transport liquid medium to preserve the viability of virus or bacteria during transport, as VTM/ UTM, and biological reagent for specimen analysis, and bacteria culture media for viewing growth.
- The transparent polypropylene tubes allow unobstructed viewing of the sample.
- The tubes are shatter-proof and may be exposed to freezing temperatures.
- The tubes are self-standing with conical bottom.
- The tubes can be supplied with a label on request
- Gamma Sterile




Swab Capture Design



| P/N | VOLUME | DIM(MM) | GRADUATIONS | TUBE/CAP | PACKING |
|------------------|--------|-----------|-------------------|----------|---|
| 1 PLAC4004X1576S | 5ml | 16 x 56.5 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 2 PLAC4004X1577S | 8ml | 16 x 85 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 3 PLAC4004X1578S | 10ml | 16 x 94 | Molded graduation | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |
| 4 PLAC4004X1593S | 12ml | 16 x 102 | Non graduated | Separate | 100tubes/bag+1,000caps/bag, 1,000tubes+1,000caps/carton |

*Dimensions for tube only




| | A(mm) | B(mm) | C(mm) | D(mm) |
|---|-------|-------|-------|-------|
| 1 | 21 | 16 | 58 | 56.5 |
| 2 | 21 | 16 | 89.5 | 85 |
| 3 | 21 | 16 | 96 | 94 |
| 4 | 21 | 16 | 102 | 100 |

Pipette Tips

Pipette Tips

Automated pipette tips are compatible with a variety of automated pipetting workstations and automated sampling systems. They are used to distribute and transfer liquids to help complete high-throughput operation of biological samples. Universal pipette tips are made with high precision molds. With excellent processing technology and good pipetting performance, they are adapted to major brands such as DragonLab, Gilson, Eppendorf, Thermofisher, etc.



- Made of high quality PP(polypropylene)
- DNase/RNase free, Pyrogen free
- Low CV Accuracy, strong hydrophobicity, no liquid adhesion
- High quality filter, tight sealing, cross infection prevention
- Excellent toughness, not easy to deform, good verticality and air tightness


Agilent



30µl

250µl

| Sort | Name | Item No. | Packing |
|-------|---|---------------|-------------------------------------|
| 30µl | AG 30µl (384 wells) transparent pipette tip | PLAKGAT030ATP | 384 pcs/box(1 rack/box), 50 box/ctn |
| | AG 30µl (384 wells) transparent pipette tip with filter | PLAKGAF030ATP | 384 pcs/box(1 rack/box), 50 box/ctn |
| 70µl | AG 70µl (384 wells) transparent pipette tip | PLAKGAT070ATP | 384 pcs/box(1 rack/box), 50 box/ctn |
| | AG 70µl (384 wells) transparent pipette tip with filter | PLAKGAF070ATP | 384 pcs/box(1 rack/box), 50 box/ctn |
| 250µl | AG 250µl (96 wells) transparent pipette tip | PLAKGAT250ATP | 96 pcs/box(1 rack/box), 50 box/ctn |
| | AG 250µl (96 wells) transparent pipette tip with filter | PLAKGAF250ATP | 96 pcs/box(1 rack/box), 50 box/ctn |



50µl

200µl

1000µl

MCA50µl

MCA200µl

| Sort | Name | Item No. | Packing |
|--|---|-------------------------------------|-------------------------------------|
| 20µl | TC 20µl transparent pipette tip | PLAKGAT020TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT020TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 20µl transparent pipette tip with filter | PLAKGAF020TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAF020TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 20µl conductive pipette tip | PLAKGAT020TB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT020TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| TC 20µl conductive pipette tip with filter | PLAKGAF020TB | 192 pcs/box(2 rack/box), 24 box/ctn | |
| | PLAKGAF020TP | 96 pcs/box (1 rack/box), 50 box/ctn | |
| 50µl | TC 50µl transparent pipette tip | PLAKGAT050TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT050TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl transparent pipette tip with filter | PLAKGAF050TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAF050TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl conductive pipette tip | PLAKGAT050TB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT050TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl conductive pipette tip with filter | PLAKGAF050TB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAF050TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl slender transparent pipette tip | PLAKGAT050TTPLP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl slender transparent pipette tip with filter | PLAKGAF050TTPLP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl slender conductive pipette tip | PLAKGAT050TLP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 50µl slender conductive pipette tip with filter | PLAKGAF050TLP | 96 pcs/box (1 rack/box), 50 box/ctn |


| Sort | Name | Item No. | Packing |
|---|--|-------------------------------------|--------------------------------------|
| 200µl | TC 200µl transparent pipette tip | PLAKGAT200TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT200TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 200µl transparent pipette tip with filter | PLAKGAF200TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAF200TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 200µl conductive pipette tip | PLAKGAT200TB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT200TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| TC 200µl conductive pipette tip with filter | PLAKGAF200TB | 192 pcs/box(2 rack/box), 24 box/ctn | |
| | PLAKGAF200TP | 96 pcs/box (1 rack/box), 50 box/ctn | |
| 1000µl | TC 1000µl transparent pipette tip | PLAKGAT1000TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT1000TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 1000µl transparent pipette tip with filter | PLAKGAF1000TTPB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAF1000TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 200µl conductive pipette tip | PLAKGAT1000TB | 192 pcs/box(2 rack/box), 24 box/ctn |
| | | PLAKGAT1000TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| TC 200µl conductive pipette tip with filter | PLAKGAF1000TB | 192 pcs/box(2 rack/box), 24 box/ctn | |
| | PLAKGAF1000TP | 96 pcs/box (1 rack/box), 50 box/ctn | |
| 5ml | TC 5ml transparent pipette tip | PLAKGAT5000TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 5ml transparent pipette tip with filter | PLAKGAF5000TTPP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 5ml conductive pipette tip | PLAKGAT5000TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC 5ml conductive pipette tip with filter | PLAKGAF5000TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| MCA15µl | TC MCA 15µl (384 wells) transparent pipette tip | PLAKGAT15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 15µl (384 wells) transparent pipette tip with filter | PLAKGAF15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| MCA50µl | TC MCA 50µl (96 wells) transparent pipette tip | PLAKGAT15M3TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 50µl (96 wells) transparent pipette tip with filter | PLAKGAF15M3TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 50µl (384 wells) transparent pipette tip | PLAKGAT15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 50µl (384 wells) transparent pipette tip with filter | PLAKGAF15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| MCA125µl | TC MCA 125µl (384 wells) transparent pipette tip | PLAKGAT15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 125µl (384 wells) transparent pipette tip with filter | CPLAKGAF15M3TP | 384 pcs/box (1 rack/box), 50 box/ctn |
| MCA200µl | TC MCA 200µl (96 wells) transparent pipette tip | PLAKGAT15M3TP | 96 pcs/box (1 rack/box), 50 box/ctn |
| | TC MCA 200µl (96 wells) transparent pipette tip with filter | PLAKGAF15M3TP | 96 pcs/box (1 rack/box), 50 box/ctn |

Hamilton



| Sort | Name | Item No. | Packing | |
|--|---|----------------------------------|---------------------------------|---------------------------------|
| 50µl | HM 50µl transparent pipette tip | PLAKGAT050HTPB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAT050HTPP | 96pcs/box(1rack/box), 50box/ctn | |
| | HM 50µl transparent pipette tip with filter | PLAKGAF050HTPB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAF050HTPP | 96pcs/box(1rack/box), 50box/ctn | |
| | HM 50µl conductive pipette tip | PLAKGAT050HB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAT050HP | 96pcs/box(1rack/box), 50box/ctn | |
| | HM 50µl conductive pipette tip with filter | PLAKGAF050HB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAF050HP | 96pcs/box(1rack/box), 50box/ctn | |
| | HM 50µl (II) conductive pipette tip | PLAKGAT0502HP | 96pcs/box(1rack/box), 50box/ctn | |
| | HM 50µl (II) conductive pipette tip with filter | PLAKGAF0502HP | 96pcs/box(1rack/box), 50box/ctn | |
| | 300µl | HM 300µl transparent pipette tip | PLAKGAT300HTPB | 480pcs/box(5rack/box), 9box/ctn |
| | | | PLAKGAT300HTPP | 96pcs/box(1rack/box), 50box/ctn |
| HM 300µl transparent pipette tip with filter | | PLAKGAF300HTPB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAF300HTPP | 96pcs/box(1rack/box), 50box/ctn | |
| HM 300µl conductive pipette tip | | PLAKGAT300HB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAT300HP | 96pcs/box(1rack/box), 50box/ctn | |
| HM 300µl conductive pipette tip with filter | | PLAKGAF300HB | 480pcs/box(5rack/box), 9box/ctn | |
| | | PLAKGAF300HP | 96pcs/box(1rack/box), 50box/ctn | |

| Sort | Name | Item No. | Packing |
|-----------------------------------|--|-----------------|---------------------------------|
| 300µl (Extended Length) | HM 300µl (Extended Length) transparent pipette tip | PLAKGAT300HTPLB | 480pcs/box(5rack/box), 9box/ctn |
| | HM 300µl (Extended Length) transparent pipette tip with filter | PLAKGAF300HTPLB | 480pcs/box(5rack/box), 9box/ctn |
| | HM 300µl (Extended Length) conductive pipette tip | PLAKGAT300HLP | 96pcs/box(1rack/box), 50box/ctn |
| | HM 300µl (Extended Length) conductive pipette tip with filter | PLAKGAF300HLP | 96pcs/box(1rack/box), 50box/ctn |
| 1000µl | HM 1000µl transparent pipette tip | PLAKGAT1000HTPB | 480pcs/box(5rack/box), 9box/ctn |
| | | PLAKGAT1000HTPP | 96pcs/box(1rack/box), 50box/ctn |
| | HM 1000µl transparent pipette tip with filter | PLAKGAF1000HTPB | 480pcs/box(5rack/box), 9box/ctn |
| | | PLAKGAF1000HTPP | 96pcs/box(1rack/box), 50box/ctn |
| | HM 1000µl conductive pipette tip | PLAKGAT1000HB | 480pcs/box(5rack/box), 9box/ctn |
| | | PLAKGAT1000HP | 96pcs/box(1rack/box), 50box/ctn |
| | HM 1000µl conductive pipette tip with filter | PLAKGAF1000HB | 480pcs/box(5rack/box), 9box/ctn |
| | | PLAKGAF1000HP | 96pcs/box(1rack/box), 50box/ctn |




| Sort | Name | Item No. | Packing |
|--------|--|----------------|---------------------------------|
| 20µl | BKM 20µl transparent pipette tip | PLAKGAT020BTP | 96pcs/box(1rack/box), 50box/ctn |
| | BKM 20µl transparent pipette tip with filter | PLAKGAF020BTP | 96pcs/box(1rack/box), 50box/ctn |
| 50µl | BKM 50µl transparent pipette tip | PLAKGAT050BTP | 96pcs/box(1rack/box), 50box/ctn |
| | BKM 50µl transparent pipette tip with filter | PLAKGAF050BTP | 96pcs/box(1rack/box), 50box/ctn |
| 250µl | BKM 250µl transparent pipette tip | PLAKGAT250BTP | 96pcs/box(1rack/box), 50box/ctn |
| | BKM 250µl transparent pipette tip with filter | PLAKGAF250BTP | 96pcs/box(1rack/box), 50box/ctn |
| 1000µl | BKM 1000µl transparent pipette tip | PLAKGAT1000BTP | 96pcs/box(1rack/box), 50box/ctn |
| | BKM 1000µl transparent pipette tip with filter | PLAKGAF1000BTP | 96pcs/box(1rack/box), 50box/ctn |

Xantus、Apricot Designs




1000µl



50µl

| Sort | Name | Item No. | Packing |
|--------|--|----------------|--------------------------------------|
| 1000µl | XTS 1000µl transparent pipette tip | PLAKGAT1000XTP | 960 pcs/box(10 rack/box), 10 box/ctn |
| | | | 96 pcs/box(4 rack/box), 20 box/ctn |
| | XTS 1000µl transparent pipette tip with filter | PLAKGAF1000XTP | 96pcs/box(1rack/box), 45box/ctn |
| 50µl | MX 50µl (96 wells) transparent pipette tip | PLAKGAT50MXTP | 96pcs/box(1rack/box), 50box/ctn |
| | MX 50µl (384 wells) transparent pipette tip | PLAKGAT50MXTP3 | 384pcs/box(1rack/box), 50box/ctn |


Universal Pipette Tip Series



| Sort | Name | Item No. | Packing |
|------------------------|---|---------------|--------------------------|
| 10µl | 10µl universal pipette tip | PLAKGPT10TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT10TP9 | 96 pcs/box, 50 box/ctn |
| | 10µl universal pipette tip with filter | PLAKGFT10TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT10TP9 | 96 pcs/box, 50 box/ctn |
| 10µl (Extended Length) | 10µl extended length (50mm) universal pipette tip | PLAKGPT10TPL | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT10TPL9 | 96 pcs/box, 50 box/ctn |
| | 10µl extended length (50mm) universal pipette tip with filter | PLAKGFT10TPL | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT10TPL9 | 96 pcs/box, 50 box/ctn |
| 20µl | 20µl universal pipette tip with filter | PLAKGFT20TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT20TP9 | 96 pcs/box, 50 box/ctn |
| 50µl | 50µl universal pipette tip with filter | PLAKGFT50TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT50TP9 | 96 pcs/box, 50 box/ctn |
| 100µl | 100µl universal pipette tip with filter | PLAKGFT100TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT100TP9 | 96 pcs/box, 50 box/ctn |
| 200µl | 200µl universal pipette tip | PLAKGPT200TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT200TP9 | 96 pcs/box, 50 box/ctn |
| | 200µl universal pipette tip with filter | PLAKGFT200TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT200TP9 | 96 pcs/box, 50 box/ctn |

| Sort | Name | Item No. | Packing |
|---|--|-------------------------|--------------------------|
| 200µl (Extended Length) | 200µl extended length (89mm) universal pipette tip | PLAKGPT200TPL | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT200TPL9 | 96 pcs/box, 50 box/ctn |
| | 200µl extended length (89mm) universal pipette tip with filter | PLAKGFT200TPL | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT200TPL9 | 96 pcs/box, 50 box/ctn |
| | 200µl extended length (89mm) Wide-mouth universal pipette tip | PLAKGPT200KTPL | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT200KTPL9 | 96 pcs/box, 50 box/ctn |
| 200µl extended length (89mm) Wide-mouth universal pipette tip with filter | PLAKGFT200KTPL | 1000 pcs/bag, 5 bag/ctn | |
| | PLAKGFT200KTPL9 | 96 pcs/box, 50 box/ctn | |
| 300µl | 300µl universal pipette tip | PLAKGFT300TP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT300TP9 | 96 pcs/box, 50 box/ctn |
| 1000µl | 1000µl universal pipette tip | PLAKGFT1000TP | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT1000TP9 | 96 pcs/box, 50 box/ctn |
| | 1000µl universal pipette tip with filter | PLAKGFT1000TP | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT1000TP9 | 96 pcs/box, 50 box/ctn |
| | 1000µl Blue universal pipette tip | PLAKGPT1000B | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT1000B9 | 96 pcs/box, 50 box/ctn |
| | 1000µl Blue universal pipette tip with filter | PLAKGFT1000B | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT1000B9 | 96 pcs/box, 50 box/ctn |
| | 1000µl Wide-mouth universal pipette tip | PLAKGPT1000KTP | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT1000KTP9 | 96 pcs/box, 50 box/ctn |
| 1000µl Wide-mouth universal pipette tip with filter | PLAKGFT1000KTP | 1000 pcs/bag, 5 bag/ctn | |
| | PLAKGFT1000KTP9 | 96 pcs/box, 50 box/ctn | |
| 1000µl (Extended Length) | 1000µl extended length (102mm) universal pipette tip | PLAKGPT1000TPL | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT1000TPL9 | 96 pcs/box, 50 box/ctn |
| | 1000µl extended length (102mm) universal pipette tip with filter | PLAKGFT1000TPL | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT1000TPL9 | 96 pcs/box, 50 box/ctn |

Rainin Pipette Tips



| Sort | Name | Item No. | Packing |
|--------|--|-----------------|--------------------------|
| 20µl | 20µl universal pipette tip | PLAKGPT20RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT20RTP9 | 96 pcs/box, 50 box/ctn |
| | 20µl universal pipette tip with filter | PLAKGFT20RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT20RTP9 | 96 pcs/box, 50 box/ctn |
| 200µl | 200µl universal pipette tip | PLAKGPT200RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT200RTP9 | 96 pcs/box, 50 box/ctn |
| | 200µl universal pipette tip with filter | PLAKGFT200RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT200RTP9 | 96 pcs/box, 50 box/ctn |
| 300µl | 300µl universal pipette tip | PLAKGPT300RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGPT300RTP9 | 96 pcs/box, 50 box/ctn |
| | 300µl universal pipette tip with filter | PLAKGFT300RTP | 1000 pcs/bag, 20 bag/ctn |
| | | PLAKGFT300RTP9 | 96 pcs/box, 50 box/ctn |
| 1000µl | 1000µl universal pipette tip | PLAKGPT1000RTP | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGPT1000RTP9 | 96 pcs/box, 50 box/ctn |
| | 1000µl universal pipette tip with filter | PLAKGFT1000RTP | 1000 pcs/bag, 5 bag/ctn |
| | | PLAKGFT1000RTP9 | 96 pcs/box, 50 box/ctn |

HLF Pipette Tips



| Sort | Name | Item No. | Packing |
|--------|---|------------------|------------------------|
| 1000µl | S-H11A 1000ul universal pipette tip with filter | PLAKGPT1000TPL9S | 96 pcs/box, 50 box/ctn |



WORLDWIDE

EUROPE

Italy Office

Headquarters

GVS S.p.A.
Via Roma 50
40069 Zola Predosa (BO) - Italy
Tel. +39 051 6176311
gvs@gvs.com

GVS SUD srl
Via Tavernole, 8, 83030 Manocalzati AV

GVS - G3
Via Guido Rossa 2/4
40069 Zona Industriale - Zola Predosa
(Bologna) Italy

Haemotronic S.p.A.
Via Carreri 16
41037 Mirandola
(Modena) Italy

Haemotronic S.p.A.
Via Ugo Roncada 83
Località Carbonara di Po

United Kingdom

GVS Filter Technology UK
Vickers Industrial Estate
Mellishaw Lane, Morecambe Lancashire
LA3 3EN
Tel. +44 (0) 1524 847600
gvsuk@gvs.com

Romania

GVS Microfiltrazione srl
Sat Ciorani de Sus 1E - Comuna Ciorani
Prahova România
Tel. (+40) 244 463044
gvsro@gvs.com

ASIA

China

GVS Technology (Suzhou) Co., Ltd.
Fengqiao Civil-Run Sci-Tech Park,
602 Changjiang Road, S.N.D.
Suzhou, China 215129
Tel. +86 512 6661 9880
gvschina@gvs.com

GVS YIBO Medical Devices Co. Ltd.
17, Zhongshan East - Yuyao city,
315403 Zhejiang Province, China
Tel. +86 574 6257 5697

Shanghai Transfusion Technology Co., Ltd.
No.500 Youdong Rd
(Shanghai), China
Tel:+86 21 54885001

Suzhou Laishi Transfusion Equipment Co., Ltd.
Changsheng Road, Tongli Town, Wujiang District
Suzhou (Jiangsu), China
Tel:+86 512-63331782
laishi@gvs.com

AMERICA

U.S.A.

GVS North America
63 Community Drive
Sanford, ME 04073 - USA
Tel. +1 866 7361250
gvsusa@gvs.com

GVS Filtration Inc.
2150 Industrial Drive
Findlay, OH. 45840 - USA
Tel. +1.419.423.9040
gvsfiltration@gvs.com

2200 W 20th Avenue
Bloomer, WI 54724 - USA
Tel. +1.715.568.5944
gvsfiltration@gvs.com

GVS Westborough
125 Flanders Rd, Westborough,
MA 01581, USA

RPB Safety LLC.
2807 Samoset Road,
Royal Oak, MI
48073 USA
tel. 1-866-494-4599

Puerto Rico

GVS Puerto Rico, LLC
98 Carr 194 - Fajardo,
Puerto Rico, 00738-2988, USA
Tel. +1.787.355.4100
gvspuertorico@gvs.com

México

GVS Filter Technology de Mexico
Universal No. 550, Vynmsa Aeropuerto Apodaca
Industrial Park, Ciudad Apodaca, Nuevo León,
C.P. 66626 - México
Tel. +52 81 2282 9003
gvs mex@gvs.com

Haemotronic de Mexico
Avenida Mike Allen, Parque Industrial Reynosa
1011-A, 88780 Reynosa, Tamaulipas., Messico

Brazil

GVS do Brasil Ltda.
Rodovia Conego Cyriaco Scaranello Pires 251
Jardim Chapadão, CEP 13193-580
Monte Mor (SP) - Brasil
Tel. +55 19 38797200
gvs@gvs.com.br

PRODUCT COLLECTION - LIFE SCIENCE PRODUCTS Plastic consumables

Copyright © 2023 GVS © S.p.A.

All Right Reserved

Printing History: Version 21062023

While every precaution has been taken in the preparation of this catalog, data are subject to change without notice. Results in specific application of GVS products may vary according to the conditions and applications. GVS assumes no responsibility for damage resulting from incorrect use of our products.