

Dispenser Tips



Dispenser Tips for Repeating Pipettor

Manufactured from high-quality virgin plastics, these tips are compatible with most standard repeating pipettes and are suitable for a wide range of applications. They are produced without chemical additives to prevent reagent contamination.

Features and Benefits

- Precision-engineered seal between piston and cylinder ensures smooth operation, enhanced accuracy, and repeatability
- Extended piston design reduces dead volume and maximizes fluid usage
- Ergonomic funnel-shaped design enables easy handling and protects gloves from damage
- Manufactured entirely from premium virgin high-density polyethylene (HDPE) and polypropylene (PP)
- Utilizes the direct displacement method, offering a safe and reliable solution for dispensing viscous, infectious, and volatile liquids
- Every tip is individually tested during automated production to ensure consistent quality and performance

Ordering information

Product Code	Description	PS/BOX
PLAKG020XTPSA	200µl Dispenser tips, sterile, 10pcs/bag,500pcs/carton	500
PLAKG01XTPXSA	1000µl Dispenser tips, sterile, 10pcs/bag,500pcs/carton	500
PLAKG5XTPXSA	5ml Dispenser tips, sterile, 10pcs/bag,500pcs/carton	500
PLAKG10XTPXSA	10ml Dispenser tips, sterile, 10pcs/bag,500pcs/carton	300
PLAKG50XTPXSA	50ml Dispenser tips, sterile, 10pcs/bag,500pcs/carton	100
PLAKG020XTPA	200µl Dispenser tips, sterile, 1pc/bag,500pcs/carton	500
PLAKG01XTPA	1000µl Dispenser tips, sterile, 1pc/bag,500pcs/carton	500
PLAKG5XTPA	5ml Dispenser tips, sterile, 1pc/bag,500pcs/carton	500
PLAKG10XTPA	10ml Dispenser tips, sterile, 1pc/bag,500pcs/carton	300
PLAKG50XTPA	50ml Dispenser tips, sterile, 1pc/bag,500pcs/carton	100

