

Pipette Stand



GVS's Scientific pipette stand holds up to six single or multichannel pipettes, protecting them from contamination. It is detachable and compatible with all Major pipette Brands.

Features

- Compact design; Occupies less bench space
- Compatible with all Major Pipette Brands
- Hold up to 6 single and multi-channel pipettes
- Special Notch Design ensures Compatibility & Easy Retrieval

Ordering information

Product Code	FLPFGPS030A
Positions	6
Color	Grey
Material	ABS
Dimension(L*W*H)	352*156*305mm
Net Weight	0.8kg

