

## Solution Bottle

GVS solution bottles are ideal for specialty chemicals, diagnostics, biologicals, reagents, adhesives and veterinary pharmaceuticals. HDPE packaging bottles are excellent for storage, shipping and packaging of liquid reagents.



- Lab quality wide-mouth and narrow-mouth bottles.
- Liner less caps provide the ultimate in leakproof protection
- Ideal for storing, shipping and packaging liquids
- HDPE/PP bottles available
- Available with 8ml/15ml/30ml/60ml/125ml/250ml/500ml/1000ml
- Natural and amber are available
- Leakage test: vacuum, 95KPA 30mins

| Product Code | Description                         | Material | Color   | Pcs/Case |
|--------------|-------------------------------------|----------|---------|----------|
| SBPGWMCH0008 | wide mouth, 8ml packaging bottle    | HDPE     | Natural | 2000     |
| SBPGWMCH0015 | wide mouth, 15ml packaging bottle   | HDPE     | Natural | 1500     |
| SBPGWMCH0030 | wide mouth, 30ml packaging bottle   | HDPE     | Natural | 800      |
| SBPGWMCH0060 | wide mouth, 60ml packaging bottle   | HDPE     | Natural | 800      |
| SBPGWMCH0125 | wide mouth, 125ml packaging bottle  | HDPE     | Natural | 500      |
| SBPGWMCH0250 | wide mouth, 250ml packaging bottle  | HDPE     | Natural | 250      |
| SBPGWMCH0500 | wide mouth, 500ml packaging bottle  | HDPE     | Natural | 125      |
| SBPGWMCH1000 | wide mouth, 1000ml packaging bottle | HDPE     | Natural | 50       |
| SBPGWMAH0008 | wide mouth, 8ml packaging bottle    | HDPE     | Amber   | 2000     |
| SBPGWMAH0015 | wide mouth, 15ml packaging bottle   | HDPE     | Amber   | 1500     |
| SBPGWMAH0030 | wide mouth, 30ml packaging bottle   | HDPE     | Amber   | 800      |
| SBPGWMAH0060 | wide mouth, 60ml packaging bottle   | HDPE     | Amber   | 800      |
| SBPGWMAH0125 | wide mouth, 125ml packaging bottle  | HDPE     | Amber   | 500      |
| SBPGWMAH0250 | wide mouth, 250ml packaging bottle  | HDPE     | Amber   | 250      |
| SBPGWMAH0500 | wide mouth, 500ml packaging bottle  | HDPE     | Amber   | 125      |
| SBPGWMAH1000 | wide mouth, 1000ml packaging bottle | HDPE     | Amber   | 50       |
| SBPGWMCP0008 | wide mouth, 8ml packaging bottle    | PP       | Natural | 2000     |
| SBPGWMCP0015 | wide mouth, 15ml packaging bottle   | PP       | Natural | 1500     |
| SBPGWMCP0030 | wide mouth, 30ml packaging bottle   | PP       | Natural | 800      |
| SBPGWMCP0060 | wide mouth, 60ml packaging bottle   | PP       | Natural | 800      |
| SBPGWMCP0125 | wide mouth, 125ml packaging bottle  | PP       | Natural | 500      |
| SBPGWMCP0250 | wide mouth, 250ml packaging bottle  | PP       | Natural | 250      |
| SBPGWMCP0500 | wide mouth, 500ml packaging bottle  | PP       | Natural | 125      |
| SBPGWMCP1000 | wide mouth, 1000ml packaging bottle | PP       | Natural | 50       |
| SBPGWMAP0008 | wide mouth, 8ml packaging bottle    | PP       | Amber   | 2000     |
| SBPGWMAP0015 | wide mouth, 15ml packaging bottle   | PP       | Amber   | 1500     |
| SBPGWMAP0030 | wide mouth, 30ml packaging bottle   | PP       | Amber   | 800      |
| SBPGWMAP0060 | wide mouth, 60ml packaging bottle   | PP       | Amber   | 800      |
| SBPGWMAP0125 | wide mouth, 125ml packaging bottle  | PP       | Amber   | 500      |
| SBPGWMAP0250 | wide mouth, 250ml packaging bottle  | PP       | Amber   | 250      |
| SBPGWMAP0500 | wide mouth, 500ml packaging bottle  | PP       | Amber   | 125      |
| SBPGWMAP1000 | wide mouth, 1000ml packaging bottle | PP       | Amber   | 50       |

| Product Code | Description                           | Material | Color   | Pcs/Case |
|--------------|---------------------------------------|----------|---------|----------|
| SBPGNMCH0030 | narrow mouth, 30ml packaging bottle   | HDPE     | Natural | 800      |
| SBPGNMCH0060 | narrow mouth, 60ml packaging bottle   | HDPE     | Natural | 800      |
| SBPGNMCH0125 | narrow mouth, 125ml packaging bottle  | HDPE     | Natural | 500      |
| SBPGNMCH0250 | narrow mouth, 250ml packaging bottle  | HDPE     | Natural | 250      |
| SBPGNMCH0500 | narrow mouth, 500ml packaging bottle  | HDPE     | Natural | 125      |
| SBPGNMCH1000 | narrow mouth, 1000ml packaging bottle | HDPE     | Natural | 50       |
| SBPGNMAH0030 | narrow mouth, 30ml packaging bottle   | HDPE     | Amber   | 800      |
| SBPGNMAH0060 | narrow mouth, 60ml packaging bottle   | HDPE     | Amber   | 800      |
| SBPGNMAH0125 | narrow mouth, 125ml packaging bottle  | HDPE     | Amber   | 500      |
| SBPGNMAH0250 | narrow mouth, 250ml packaging bottle  | HDPE     | Amber   | 250      |
| SBPGNMAH0500 | narrow mouth, 500ml packaging bottle  | HDPE     | Amber   | 125      |
| SBPGNMAH1000 | narrow mouth, 1000ml packaging bottle | HDPE     | Amber   | 50       |
| SBPGNMCP0030 | narrow mouth, 30ml packaging bottle   | PP       | Natural | 800      |
| SBPGNMCP0060 | narrow mouth, 60ml packaging bottle   | PP       | Natural | 800      |
| SBPGNMCP0125 | narrow mouth, 125ml packaging bottle  | PP       | Natural | 500      |
| SBPGNMCP0250 | narrow mouth, 250ml packaging bottle  | PP       | Natural | 250      |
| SBPGNMCP0500 | narrow mouth, 500ml packaging bottle  | PP       | Natural | 125      |
| SBPGNMCP1000 | narrow mouth, 1000ml packaging bottle | PP       | Natural | 50       |
| SBPGNMAP0030 | narrow mouth, 30ml packaging bottle   | PP       | Amber   | 800      |
| SBPGNMAP0060 | narrow mouth, 60ml packaging bottle   | PP       | Amber   | 800      |
| SBPGNMAP0125 | narrow mouth, 125ml packaging bottle  | PP       | Amber   | 500      |
| SBPGNMAP0250 | narrow mouth, 250ml packaging bottle  | PP       | Amber   | 250      |
| SBPGNMAP0500 | narrow mouth, 500ml packaging bottle  | PP       | Amber   | 125      |
| SBPGNMAP1000 | narrow mouth, 1000ml packaging bottle | PP       | Amber   | 50       |

