



Pipette Tips

Pipette Tips

GVS provides a variety on pipette tips from standard to extra-long design, in different packs of them, bulk/rack/refill, with or without filter and low retention types.

All sterile products are irradiated electron beam radiation, which is the same level as gamma irradiation but less harsh on plastics than radiation and does not leave toxic residues like gas sterilization usually does to the product.

All our sterile pipette tips are manufactured out of the highest-grade resins and are certified to be RNase, DNase, DNA, pyrogen, and ATP free. Our non-sterile tips are made from the same high-quality resins and are certified RNase and DNase free as well.



Liquid Handling

10uL, DNase/Rnase free, 32mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10RSA	0.1-10ul	-	Yes	rack pack	96/rack, 50racks/case
PLAPGF10SA	0.1-10ul	Yes	Yes	rack pack	96/rack, 50racks/case
PLAPGREF10SA	0.1-10ul	-	Yes	refill pack	96/rack, 10racks/unit, 10units/case
PLAPG10BA	0.1-10ul	-	-	bulk pack	1,000/bag, 40bags/case
PLAPGF10BA	0.1-10ul	Yes	-	bulk pack	1,000/bag, 40bags/case

10uL extra long design, DNase/Rnase free, 46mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10ERSA	0.1-20ul	-	Yes	rack pack	96/rack, 50racks/case
PLAPGF10ESA	0.1-20ul	Yes	Yes	rack pack	96/rack, 50racks/case
PLAPG10EBA	0.1-20ul	-	-	bulk pack	1,000/bag, 40bags/case
PLAPGF10EBA	0.1-20ul	Yes	-	bulk pack	1,000/bag, 40bags/case

10uL extra long, thin and sharp design, DNase/Rnase free, 46mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10ETRSA	0.1-20ul	-	Yes	rack pack	96/rack, 50racks/case
PLAPGF10ETSA	0.1-20ul	Yes	Yes	rack pack	96/rack, 50racks/case
PLAPGREF10ETSA	0.1-20ul	-	Yes	refill pack	96/rack, 9racks/unit, 10units/case
PLAPG10ETBA	0.1-20ul	-	-	bulk pack	1,000/bag, 40bags/case
PLAPGF10ETBA	0.1-20ul	Yes	-	bulk pack	1,000/bag, 40bags/case

200uL standard, bevelled, DNase/Rnase free, yellow, 51mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG200RSA	1-200ul	-	Yes	rack pack	96/rack, 50racks/case
PLAPGF20SA	1-20ul	Yes	Yes	rack pack	96/rack, 50racks/case
PLAPGREF200SA	1-200ul	-	Yes	refill pack	96/rack, 12racks/unit, 10units/case
PLAPG200BA	1-200ul	-	-	bulk pack	1,000/bag, 20bags/case
PLAPGF20BA	1-20ul	Yes	-	bulk pack	1,000/bag, 20bags/case

Liquid Handling

200uL standard wide, Bevelled, DNase/RNase free, 51mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPGF50SA	1-50ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF100SA	1-100ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF200SA	1-200ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF50BA	1-50ul	Yes	-	bulk pack	1,000/bag,20bags/case
PLAPGF100BA	1-100ul	Yes	-	bulk pack	1,000bag,20bags/case
PLAPGF200BA	1-200ul	Yes	-	bulk pack	1,000bag,20bags/case

200uL extra long or 300uL, straight, DNase/RNase free, 60mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG200ERSA	1-300ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF100ESA	1-150ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF200ESA	1-300ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGREF200ESA	1-300ul	-	Yes	refill pack	96/rack,10racks/unit,10units/case
PLAPG200EBA	1-300ul	-	-	bulk pack	1,000/bag,20bags/case
PLAPGF200EBA	1-300ul	Yes	-	bulk pack	1,000bag,20bags/case

1000uL standard, DNase/RNase free, blue, 78mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG1000RSA	100-1,000ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF1000SA	100-1,000ul	Yes	Yes	rack pack	96/rack, 50racks/case
PLAPG1000BA	100-1,000ul	-	-	bulk pack	1,000/bag,10bags/case
PLAPGF1000BA	100-1,000ul	Yes	-	bulk pack	1,000/bag,10bags/case

1000uL extra long or 1250uL, DNase/RNase free, blue, 102mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPGF1000ESA	100-1,000ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGREF1000ESA	100-1,200ul	-	Yes	refill pack	96/rack,8racks/unit,10units/case

Liquid Handling

Medium long 1000µL DNase/RNase free, blue, 88mm(Some Special Application)



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPGF1000MSA	100-1,000ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPG1000MBA	100-1,000ul	-	-	bulk pack	1,000/bag,10bags/case
PLAPGF1000MBA	100-1,000ul	Yes	-	bulk pack	1,000/bag,10bags/case

Low Retention Pipette Tips

Low Retention pipette tips are specifically designed for applications requiring high accuracy and reproducibility, especially pipetting viscous liquids in complex reactions such as PCR, RT PCR, and DNA sequencing. Pipetting viscous liquids also creates pipetting problems due to the inability to expel fluid because of surface adhesion between the liquid and the inside of the pipette tip. This fluid retention can be quite a problem because residual liquid typically remains in the pipette tip after the operation.

Features

- Improves Sample Accuracy.
- Increases Reliability.
- Saves Expensive Reagents.
- Precise Sample Delivery.
- Sterile Product is certified RNase, ATP, DNA, DNase and pyrogen free.

10uL, DNase/RNase free, low retention, 32mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10RLSA	0.1-10ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF10LSA	0.1-10ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRF10LSA	0.1-10ul	-	Yes	refill pack	96/rack,10racks/unit,10units/case
PLAPG10BLA	0.1-10ul	-	-	bulk pack	1,000/bag,40bags/case
PLAPGF10BLA	0.1-10ul	Yes	-	bulk pack	1,000/bag,40bags/case

10uL extra long design, DNase/RNase free, low retention, 46mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10ERLSA	0.1-20ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF10ELSA	0.1-20ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPG10EBLA	0.1-20ul	-	-	bulk pack	1,000/bag,40bags/case
PLAPGF10EBLA	0.1-20ul	Yes	-	bulk pack	1,000/bag,40bags/case

Liquid Handling

10uL extra long, thin and sharp design, DNase/RNase free, low retention, 46mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG10ETRLSA	0.1-20ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF10ETLSA	0.1-20ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRF10ETLSA	0.1-20ul	-	Yes	refill pack	96/rack,9racks/unit,10units/case
PLAPG10ETBLA	0.1-20ul	-	-	bulk pack	1,000/bag,40bags/case
PLAPGF10ETBLA	0.1-20ul	Yes	-	bulk pack	1,000/bag,40bags/case

200uL standard wide, bevelled, DNase/RNase free, low retention, clear, 51mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG200RLSA	1-200ul	-	-	rack pack	96/rack,50racks/case
PLAPGF20LSA	1-20ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRF200LSA	1-200ul	-	Yes	refill pack	96/rack,12racks/unit,10units/case
PLAPG200BLA	1-200ul	-	-	bulk pack	1,000/bag,20bags/case
PLAPGF20BLA	1-20ul	Yes	-	bulk pack	1,000/bag,20bags/case

200uL standard, bevelled, DNase/RNase free, low retention, 51 mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPGF50LSA	1-50ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF100LSA	1-100ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF200LSA	1-200ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF50BLA	1-50ul	Yes	-	bulk pack	1,000/bag,20bags/case
PLAPGF100BLA	1-100ul	Yes	-	bulk pack	1,000/bag,20bags/case
PLAPGF200BLA	1-200ul	Yes	-	bulk pack	1,000/bag,20bags/case

200uL extra long or 300uL, straight, DNase/RNase free, low retention, 60mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG200ERLSA	1-300ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF100ELSA	1-150ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGF200ELSA	1-300ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRF200ELSA	1-300ul	-	Yes	refill pack	96/rack,10racks/unit,10units/case
PLAPG200EBLA	1-300ul	-	-	bulk pack	1,000/bag,20bags/case
PLAPGF200EBLA	1-300ul	Yes	-	bulk pack	1,000/bag,20bags/case

Liquid Handling

1000uL standard, DNase/RNase free, low retention, 78mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG1000RLSA	100-1,000 ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF1000LSA	100-1,000ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPG1000BLA	100-1,000ul	-	-	bulk pack	1,000/bag,10bags/case
PLAPGF1000BLA	100-1,000ul	Yes	-	bulk pack	1,000/bag,10bags/case

1000uL extra long or 1250uL, DNase/RNase free, low retention, 102mm




Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG1000ERLSA	100-1,200ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF1000ELSA	100-1,000ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRF1000ELSA	100-1,200ul	-	Yes	refill pack	96/rack,8racks/unit,10units/case
PLAPG1000EBLA	100-1,200ul	-	-	bulk pack	1,000/bag,10bags/case
PLAPGF1000EBLA	100-1,000ul	Yes	-	bulk pack	1,000/bag,10bags/case

Medium long 1000uL, DNase/RNase free, low retention, 88mm



Product Code	Volume	With Filter	Sterile	Package Type	Pack
PLAPG1000MRLSA	100-1,000 ul	-	Yes	rack pack	96/rack,50racks/case
PLAPGF1000MLSA	100-1,000 ul	Yes	Yes	rack pack	96/rack,50racks/case
PLAPG1000MBLA	100-1,000ul	-	-	bulk pack	1,000/bag,10bags/case
PLAPGF1000MBLA	100-1,000ul	Yes	-	bulk pack	1,000/bag,10bags/case



Pipette Tip (To Fit LTS Pipette)

Features

- Easy to handle for operator
- Easier/lighter to eject
- More accurate sample delivery.
- Seal well between Pipette and Pipette tip
- Available in rack/refill pack with or without filter
- Low retention version is available

20uL Fit LTS Pipette, sterile



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT20RSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT20RLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT20RFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT20RFLSA	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRT20REFSA	-	-	refill pack	96/rack,10racks/unit,10units/case
PLAPGRT20REFLSA	-	Yes	refill pack	96/rack,10racks/unit,10units/case

200uL Fit LTS Pipette, sterile



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT200RSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT200RLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT200RFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT200RFLSA	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRT200BA	-	-	bulk pack	1,000/bag,20bags/case
PLAPGRT200BLA	-	Yes	bulk pack	1,000/bag,20bags/case
PLAPGRT200BFA	Yes	-	bulk pack	1,000/bag,20bags/case
PLAPGRT200BFLA	Yes	Yes	bulk pack	1,000/bag,20bags/case
PLAPGRT200REFSA	-	-	refill pack	96/rack,10racks/unit,10units/case
PLAPGRT200REFLSA	-	Yes	refill pack	96/rack,10racks/unit,10units/case

300uL Fit LTS Pipette, sterile



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT300RSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT300RLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT300RFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT300RFLSA	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRT300BA	-	-	bulk pack	1,000/bag,20bags/case
PLAPGRT300BLA	-	Yes	bulk pack	1,000/bag,20bags/case
PLAPGRT300BFA	Yes	-	bulk pack	1,000/bag,20bags/case
PLAPGRT300BFLA	Yes	Yes	bulk pack	1,000/bag,20bags/case
PLAPGRT300REFSA	-	-	refill pack	96/rack,8racks/unit,10units/case
PLAPGRT300REFLSA	-	Yes	refill pack	96/rack,8racks/unit,10units/case

1000uL Fit LTS Pipette, sterile




Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1000RSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT1000RLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT1000RFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT1000RFLSA	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRT1000BA	-	-	bulk pack	1,000/bag,10bags/case
PLAPGRT1000BLA	-	Yes	bulk pack	1,000/bag,10bags/case
PLAPGRT1000BFA	Yes	-	bulk pack	1,000/bag,10bags/case
PLAPGRT1000BFLA	Yes	Yes	bulk pack	1,000/bag,10bags/case
PLAPGRT1000REFSA	-	-	refill pack	96/rack,8racks/unit,10units/case
PLAPGRT1000REFLSA	-	Yes	refill pack	96/rack,8racks/unit,10units/case

1200uL Fit LTS Pipette, sterile



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1200RSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT1200RLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT1200RFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT1200RFLSA	Yes	Yes	rack pack	96/rack,50racks/case
PLAPGRT1200BA	-	-	bulk pack	1,000/bag,10bags/case
PLAPGRT1200BLA	-	Yes	bulk pack	1,000/bag,10bags/case
PLAPGRT1200BFA	Yes	-	bulk pack	1,000/bag,10bags/case
PLAPGRT1200BFLA	Yes	Yes	bulk pack	1,000/bag,10bags/case



Pipette Tip (Large Volume Tips)

- Easy to Use: Soft and flexible material for easy tip loading & ejection, universal filtered pipette tips lower the risk of repetitive stress injury (RSI).
- The smooth inner surface can prevent leakage and residues.
- Convenient & Sterile: Pipette tips rack pack allows direct loading into pipettor without hand touch, higher sterility and efficiency as well as a relief for repeated acts.
- Widely Application Compatible with most brands such as Eppendorf, Biohit, DLAB, Brand, Thermo, etc.

5ml, clear, fit Eppendorf, Biohit, DLAB Pipette



Product Code	With Filter	Sterile	Package Type	Pack
PLAPGPT5000RESA	-	Yes	rack pack	24/rack, 30 racks/case
PLAPGPT5000RFESA	Yes	Yes	rack pack	24/rack, 30 racks/case
PLAPGPT5000BEA	-	-	bulk pack	100/bag, 20 bags/case
PLAPGPT5000BFEA	Yes	-	bulk pack	100/bag, 20 bags/case

10ml, clear, fit Eppendorf, Biohit, DLAB Pipette



Product Code	With Filter	Sterile	Package Type	Pack
PLAPGPT10000RESA	-	Yes	rack pack	24/rack, 40 racks/case
PLAPGPT10000RFESA	Yes	Yes	rack pack	24/rack, 40 racks/case
PLAPGPT10000BEA	-	-	bulk pack	100/bag, 10 bags/case
PLAPGPT10000BFEA	Yes	-	bulk pack	100/bag, 10 bags/case

5ml, clear, fit Thermo, Brand Pipette



Product Code	With Filter	Sterile	Package Type	Pack
PLAPGPT5000RTSA	-	Yes	rack pack	40/rack, 40 racks/case
PLAPGPT5000RFTSA	Yes	Yes	rack pack	40/rack, 40 racks/case
PLAPGPT5000BTA	-	-	bulk pack	100/bag, 20 bags/case
PLAPGPT5000BFTA	Yes	-	bulk pack	100/bag, 20 bags/case

10ml, clear, fit Thermo, Brand Pipette



Product Code	With Filter	Sterile	Package Type	Pack
PLAPGPT10000RTSA	-	Yes	rack pack	24/rack, 40 racks/case
PLAPGPT10000RFTSA	Yes	Yes	rack pack	24/rack, 40 racks/case
PLAPGPT10000BTA	-	-	bulk pack	100/bag, 10 bags/case
PLAPGPT10000BFTA	Yes	-	bulk pack	100/bag, 10 bags/case

Robotic Pipette Tip

GVS offers a wide range of pipet tips for automated robotic workstations like, Tecan, Hamilton, Beckman etc.

All our sterile pipette tips are manufactured out of the highest-grade resins and are certified to be RNase, DNase, DNA, pyrogen, and ATP free. Our non-sterile tips are made from the same high-quality resins and are certified RNase and DNase free as well.

Automated workstations require much tighter tolerances to robotic tips than those of standard manual pipet tips. So, we strongly suggest purchasing pre-sterilized products than non-sterilized products which needs autoclaving, and autoclaving can warp the tips slightly.

50uL Clear Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT050TSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT050TLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT050TFSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT050TFLSA	Yes	Yes	rack pack	96/rack,30racks/case

200uL Clear Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT200TSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT200TLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT200TFSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT200TFLSA	Yes	Yes	rack pack	96/rack,30racks/case

1000uL Clear Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1000TSA	-	-	rack pack	96/rack,60racks/case
PLAPGRT1000TLSA	-	Yes	rack pack	96/rack,60racks/case
PLAPGRT1000TFSA	Yes	-	rack pack	96/rack,60racks/case
PLAPGRT1000TFLSA	Yes	Yes	rack pack	96/rack,60racks/case

50uL Conductive Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT050TXSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT050TXLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT050TFXSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT050TFXLSA	Yes	Yes	rack pack	96/rack,30racks/case

200uL Conductive Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT200TXSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT200TXLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT200TFXSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT200TFXLSA	Yes	Yes	rack pack	96/rack,30racks/case

1000uL Conductive Tip, sterile, Fit Robotic Tip Tecan



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1000TXSA	-	-	rack pack	96/rack,60racks/case
PLAPGRT1000TXLSA	-	Yes	rack pack	96/rack,60racks/case
PLAPGRT1000TFXSA	Yes	-	rack pack	96/rack,60racks/case
PLAPGRT1000TFXLSA	Yes	Yes	rack pack	96/rack,60racks/case

50uL Clear Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT050HSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT050HLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT050HFSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT050HFLSA	Yes	Yes	rack pack	96/rack,30racks/case

300uL Clear Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT300HSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT300HLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT300HFSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT300HFLSA	Yes	Yes	rack pack	96/rack,30racks/case

1000uL Clear Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1000HSA	-	-	rack pack	96/rack,60racks/case
PLAPGRT1000HLSA	-	Yes	rack pack	96/rack,60racks/case
PLAPGRT1000HFSA	Yes	-	rack pack	96/rack,60racks/case
PLAPGRT1000HFLSA	Yes	Yes	rack pack	96/rack,60racks/case

50uL Conductive Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT050HXSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT050HXLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT050HFXSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT050HFXLSA	Yes	Yes	rack pack	96/rack,30racks/case

300uL Conductive Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT300HXSA	-	-	rack pack	96/rack,30racks/case
PLAPGRT300HXLSA	-	Yes	rack pack	96/rack,30racks/case
PLAPGRT300HFXSA	Yes	-	rack pack	96/rack,30racks/case
PLAPGRT300HFXLSA	Yes	Yes	rack pack	96/rack,30racks/case

1000uL Conductive Tip, sterile, Fit Robotic Tip Hamilton



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT1000HXSA	-	-	rack pack	96/rack,60racks/case
PLAPGRT1000HXLSA	-	Yes	rack pack	96/rack,60racks/case
PLAPGRT1000HFXSA	Yes	-	rack pack	96/rack,60racks/case
PLAPGRT1000HFXLSA	Yes	Yes	rack pack	96/rack,60racks/case

50uL Clear Tip, rack, sterile, Fit Robotic Tip Beckman



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT050BSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT050BLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT050BFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT050BFLSA	Yes	Yes	rack pack	96/rack,50racks/case

250uL Clear Tip, rack, sterile, Fit Robotic Tip Beckman



Product Code	With Filter	Low Retention	Package Type	Pack
PLAPGRT250BSA	-	-	rack pack	96/rack,50racks/case
PLAPGRT250BLSA	-	Yes	rack pack	96/rack,50racks/case
PLAPGRT250BFSA	Yes	-	rack pack	96/rack,50racks/case
PLAPGRT250BFLSA	Yes	Yes	rack pack	96/rack,50racks/case

Note:

Fit Robotic Tip Tecan: Fits Tecan Genesis Freedom, Freedom EVO and Miniprep with LiHa

Fit Robotic Tip Hamilton: Fits Hamilton Star, Starlet and Nimbus

Fit Robotic Tip Beckman: Fits Beckman FX/NX, Multimek AP96 and Biomek3000

Robotic Tip


These premium pipette tips deliver precise performance for automated systems:

- Broad Compatibility: Designed for Fontus® and JANUs® G3 VariSpan™ modules
- Conductive Material: Graphite-infused polypropylene enables current flow




Features

- Material: PP
- Filter: Optional
- Format: Hanging
- Instrument Compatibility: Fontus Workstation and JANUS G3 Workstation
- Sterility: SAL=10⁻⁶
- Tip Type: Conductive



Ordering information

Product Code	Description	Package
PLAPGWB1000PXSA	1 mL Conductive robotic tip, wide-bore, sterile	96 pcs/rack, 24 racks/case



FLAMETIP Pipette Micro Filter Tips



Description

Natural, with filter, sterilized in rack

Purpose

Use for liquid transfer, zero retention filter.

Materials

Tips: PP (Polypropylene)

Color: Natural

Rack box: PP (Polypropylene)

Features

- FLAMETIP are manufactured from super clear high quality Polypropylene
- The surfaces of the tips are produced through a special process, this process makes the tip inner surface become super hydrophobic
- Preferred accessories for most brand micropipettor
- Tips with PP filter are also available
- Packaged in re-sealable plastic bags or extra-rigid autoclavable racks
- Every rack or case is printed with lot No. for quality traceability
- Non-Pyrogenic and validated per FDA guidelines on LAL testing for medical devices and company guidelines. The acceptance level for product is less than 0.5 EU/ml
- DNase/RNase-free
- Sterilized by gamma irradiation SAL 10-6(1S011137)
- Shelf life: 3 years after month of production
- Manufactured in a Class 100,000 cleanroom environment
- Manufactured under ISO13485:2016 and ISO9001:2015 quality management system
- This product has been tested and is free of any detectable Nucleic Acid.
- All material are certified metal-free which compliance with US Pharmacopoeia guidelines.
- This product has been tested and is free of ATP

Ordering information

Product Code	Description	Sterile	Qty/PK	Qty/Box	Qty/Case
FTS252010	10 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS233010	10 µl , long, natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252020	20 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252100	100 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS231200	200 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS233300	300 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS252000	1000 µl , natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920
FTS371000	1000 µl , long, natural, with filter, sterilized, packaged in rack box	Yes	96	960	1920